

RS-DSC



Site Preparation Guide

Ideal Setup

Table of Contents	2
Ideal Setup	3
System Components.....	4
Instrument Measurements	5
Utility Requirements.....	6
Computer Requirements	7
Accessories	8
Site Preparation Checklist	9
TA Instrument Offices.....	10



Circulator



Power



Cooling



Gas



LN₂



Fluid



Light



Hardware



Software



Temp



Lab




Customer

Ideal Setup



RS-DSC: IDEAL PLACEMENT AND BENCH MEASUREMENTS

Select a location with adequate floor space and a rigid laboratory bench that is level and is in a vibration-free environment.



Bench width: 2.4 m (8 ft)

Bench depth: 76 cm (30 in) min.

Distance from the wall: 20 cm (8 in) min.




Allow 38 cm (15 in) floor space and 20 cm (8 in) clearance at the rear of the water circulator for air circulation.

System Components



MAIN SYSTEM COMPONENTS



- A. Computer (Controller)
- B. Instrument
- C. External Water Circulator/Controller

Instrument Measurements



RS-DSC



Height LID OPEN: 65 cm (25.5 in)

Height LID CLOSED: 35.5 cm (14 in)

Width: 84 cm (33 in)

Depth: 66 cm (26 in)

Weight: 54.4 kg (120 lbs)

Utility Requirements




LABORATORY

Item	Requirement
Temperature	15–30°C
Relative Humidity	5–80% (non-condensing)
Temperature Stability	±1°C with changes in temperature being gradual shifts instead of fast changes
Instrument Location Environment	<ul style="list-style-type: none"> • Dust-free • Vibration-free • Away from exposure to direct sunlight and direct air drafts • Pollution Degree 2 Environment • Maximum altitude: 2000 m (6560 ft)



POWER

Item	Requirement
Instrument	100–240 VAC, 3.15A, 50/60 Hz
Power line	<p>Grounded, single-phase line for instrument and computer, not shared with motors, heaters, or compressors</p> <ul style="list-style-type: none"> • 15 A for voltages near 120 VAC • 10 A for voltages near 230 VAC
Electrical power cord	<ul style="list-style-type: none"> • The plug of the cord must be rated to carry at least 125% of the product's current rating. • The cord length must be less than 4.5 meters and must be UL or CSA approved.
Customer-supplied	Surge suppressor power strip



WATER

Instrument	Requirement
RS-DSC	Requires an external water circulator
TA water bath	If using a water bath from TA Instruments: 7 L (1.85 gal) distilled water

Computer Requirements



HARDWARE REQUIREMENTS

Item	Requirement
Processor	Intel® Core™ 2 Duo (2.93 GHz, 3 MB L2 cache) or better
Memory	Minimum: ≥ 16 GB RAM Recommended: ≥ 32 GB RAM
Hard drive	Minimum: ≥ 500 GB free space
CD-ROM or DVD (Optional)	DVD-ROM drive recommended if internet access is unavailable
Screen resolution	Minimum: 1024 x 768 with ≥ 24-bit colors
Graphic memory	128 MB
USB Port	Required: One available USB port (2.0 compliant)



SOFTWARE REQUIREMENTS

Item	Requirement
Operating System	<ul style="list-style-type: none">• Windows 10 or 11• Home version not supported• 64-bit version required
Browser	Internet Explorer
Service Pack	Microsoft Operating System Service Pack
Updates	Windows Operating System and associated Microsoft updates must be up to date
Network	<i>TA Instruments is not responsible for resolving issues associated with connections to your corporate network.</i>
Conflicts	<i>TA Instruments is not responsible for resolving hardware/software conflicts created by the addition of third-party hardware or software to the computer.</i>

Accessories



WATER CIRCULATOR



Height: 61.7 cm (24.3 in)

Width: 22.1 cm (8.7 in)

Depth: 54.1 cm (21.3 in)

Weight: 40.8 kg (90 lbs)

Requirements



Must be placed 20 cm (8 in) from the rear wall.



If purchased from TA, 7L (2 gal) distilled water is required for installation.



Circulator



Power



Cooling



Gas



LN₂



Fluid



Light



Hardware



Software



Temp



Lab








Customer

RS-DSC

Site Preparation Checklist



RS-DSC

	<p>Sufficient bench space for instrument, computer, sample prep space and tools:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Depth: 76 cm (30 in) minimum <input type="checkbox"/> Distance from the wall: 20 cm (8 in) minimum <p>Laboratory conditions meet the following requirements:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Dust-free, vibration-free, away from exposure to direct sunlight/air drafts, in a pollution degree 2 environment <input type="checkbox"/> Maximum altitude is 2000 m (6560 ft) <input type="checkbox"/> Relative humidity is 5–80% non-condensing <input type="checkbox"/> Temperature is 15–30°C with a stability of ±1°C
	<ul style="list-style-type: none"> <input type="checkbox"/> RS-DSC: 100–240 VAC, 4A, 50/60 Hz <input type="checkbox"/> Water Circulator: 120 VAC, 60 Hz, 13A or 240 VAC, 50 Hz, 13A
	<ul style="list-style-type: none"> <input type="checkbox"/> External water circulator <ul style="list-style-type: none"> <input type="checkbox"/> If purchased from TA Instruments, 7L (2 gal) distilled water required for installation
	<ul style="list-style-type: none"> <input type="checkbox"/> The Customer must be available to assist the TA Service Representative with lifting the instrument during installation. <input type="checkbox"/> The Customer assumes responsibility for any damage that occurs when the instrument is moved by someone other than a trained TA Instruments Service Representative.
	<ul style="list-style-type: none"> <input type="checkbox"/> Computer meets all hardware requirements <input type="checkbox"/> Computer meets all software requirements <input type="checkbox"/> Customer's IT personnel has provided Administrative privileges on the controller computer <input type="checkbox"/> The Customer's IT personnel will be on site the day of installation

I hereby acknowledge that all utility requirements have been met per the checklist above and that they will be ready at the agreed time of installation.

If all utility requirements are not met at the agreed time of installation, additional charges may be incurred for a return Service trip.

Customer _____ DD / MM / YYYY

Company _____ City _____ State _____ Country _____

Please send a signed copy of the completed checklist to your local Service representative.

TA Instruments Offices

For information on our latest products, contact information, and more, see our website at:
<http://www.tainstruments.com>.

To find your local TA Instruments office and contact information, visit
<http://www.tainstruments.com/contact/ta-directory/>

TA Instruments – Waters LLC
Corporate Headquarters
159 Lukens Drive
New Castle, DE 19720
USA

Telephone: 302-427-4000
Fax: 302-427-4001
Email: info@tainstruments.com