



TGA & SDT Theory and Applications Online Courses

Part 1: Theory & Instrumentation



©2020 TA Instruments

1

Agenda

- Understanding TGA and SDT techniques
 - Theory
 - Instrumentation
- Calibration & Verification
- Instrument & Method Considerations
 - Purge gas
 - Sample Pans
 - Sample Preparation
 - Maintenance
 - Experimental Setup
 - Experimental Methods

2



2

What is Thermogravimetric Analysis (TGA)?

- TGA measures weight/mass change (loss or gain) and the rate of weight change as a function of temperature, time and atmosphere.



3



3

What TGA Can Tell You

- Thermal Stability of Materials
- Oxidative Stability of Materials
- Composition of Multi-component Systems
- Estimated Lifetime of a Product
- Decomposition Kinetics of Materials
- The Effect of Reactive or Corrosive Atmospheres on Materials
- Moisture and Volatiles Content of Materials
- Residue

4



4

Mechanisms of Weight Change in TGA

- Weight Loss:
 - Decomposition: The breaking apart of chemical bonds.
 - Evaporation: The loss of volatiles with elevated temperature.
 - Reduction: Interaction of sample to a reducing atmosphere (hydrogen, ammonia, etc.).
 - Desorption.

- Weight Gain:
 - Oxidation: Interaction of the sample with an oxidizing atmosphere.
 - Absorption.

- All of these are kinetic processes (i.e. there is a rate at which they occur).

5



5

DSC-TGA (SDT): The Technique

- Simultaneous DSC-TGA measures both heat flow and weight changes in a material as a function of temperature or time in a controlled atmosphere from room temperature to 1500°C.
- Information obtained allows differentiation between endothermic and exothermic events which have no associated weight loss (e.g., melting and crystallization), and those which involve a weight loss (e.g., degradation).



6



6


What Simultaneous DSC-TGA Can Tell You

- Thermal Stability of Materials
- Oxidative Stability of Materials
- Composition of Multi-component Systems
- Estimated Lifetime of a Product
- Decomposition Kinetics of Materials
- The Effect of Reactive or Corrosive Atmospheres on Materials
- Moisture and Volatiles Content of Materials
- Residue
- Transition Temperatures
- Heats of Fusion and Reactions
- Melting and Boiling Points
- 7 ▪ Heat capacity



7

Thermal Stability of Polymers




8



8

Oxidative Stability (Polypropylene)




9



9

TGA of Drug A Monohydrate



Sample: Drug A Monohydrate
 Size: 15.1740mg Heating
 Rate: 10°C/min

10



10


TA Instruments TGA/SDT Models



TGA 55
TGA 550
TGA 5500



Discovery TGA



Q5000 IR



SDT 650



Q600

11



11

TGA Specifications

	TGA 5500	TGA 550/50
Temperature Range	Ambient to 1200°C	Ambient to 1000°C
Heating Rate Range	0.1 to 500°C/min (Linear) >1600°C/min (Ballistic)	0.1 to 100°C/min (Linear)
Sample Weight Capacity	1000 mg	1000 mg
Dynamic Weighing Range	1000 mg	1000 mg
Baseline Dynamic Drift (50-1000°C)	< 10 µg	<50 µg

12



12

SDT Specifications


	SDT 650
Temperature Range	Ambient to 1500°C
Heating Rate Range	0.1 to 100°C/min (Linear)
Sample Weight Capacity	200 mg
Baseline Dynamic Drift (50-1000°C) (1000°C to 1500°C)	<50 µg <50 µg

13



13

TGA: Schematic Diagram




14



14

TGA Balance and Operation



- Null-balance principle operation
- Current is applied to the meter movement
- Amount of current applied is proportional to the weight change

15



15

TGA Furnace Options: Wire Wound Furnace

- Standard furnace for TGA 55 and 550
- Ambient to 1000 °C
- Linear controlled heating rates of 0.01 to 100 °C/min
- Ballistic heating rates >600 °C/min
- Exchangeable with EGA furnace



Wire Wound (Pt/Rh) Furnace

Flow rate Balance/Sample :
40/60 ml/min

16



16

TGA Furnace Options: EGA Furnace

- Optional for TGA55 and 550
- Ambient to 1000 °C
- Linear controlled heating rates of 0.01 to 50 °C/min
- Quartz liner makes furnace easy to clean
- Exchangeable with wire wound furnace



EGA Furnace

Flow rate Balance/Sample :
10/90 ml/min

17



17

TGA Furnace Options: Infra Red (IR) Furnace

- Quartz halogen lamps as heating source
- Ambient to 1200°C
- Linear controlled heating rates of 0.01 to 500 °C/min
- Ballistic heating rates >1500 °C/min
- Integrated electromagnet for Temperature calibration with Curie point standards
- Evolved Gas Analysis capacity



IR Furnace


Flow rate Balance/Sample :
25/25 ml/min

18



18

SDT: Schematic Diagram



19



19

Calibration & Verification



20

TGA Calibration and Verification

- Two types of calibration are needed:
 - Weight
 - Temperature
- Temperature calibration is affected by:
 - Purge gas and flow rate
 - Thermal conductivity of helium \neq Thermal conductivity of nitrogen/air/oxygen \neq Thermal conductivity of argon
 - Pan type
 - Heating rates

21



21

General Calibration and Verification Guidelines

- Calibration
 - Use Calibration Mode
 - Calibrate upon installation
 - Re-calibrate if does not pass verification or if instrument setup is modified (see previous slide)
- Verification
 - Determine how often to verify data
 - Run a reference material as a sample (in standard mode)
 - Compare results vs literature values
 - Re-calibrate if results are out of tolerance

22



22

Requirements Prior to Calibration

- The TGA pan should be cleaned prior to calibration procedures.
- The purge gas flow rate should be set (see flow rates according to furnace type). The flow rate should not deviate by more than +/- 5ml/min.
- Use high purity reference materials (>99.99%) for calibration

23



23

Mass/Weight Calibration

- Weight calibration can be performed:
 - Manual, using an empty pan and calibration weights
 - Automatically, using the three weight calibration fixtures (pans) P/N 957341.901 (on Q5000IR, Discovery, and TGA 55XX).
- The calibration pans may only be used when a platinum reference pan is installed.

24



24

ASTM E 2040 - Mass Scale Calibration of Thermogravimetric Analyzers

- The mass signal generated by a TGA is compared to the mass of a reference material traceable to a national reference laboratory. A linear correlation using two calibration points is used to relate the mass (or weight) signal generated by the TGA and that of the reference material
- This test method calibrates or demonstrates conformity of thermogravimetric apparatus at ambient conditions. Most TGA experiments are carried out under temperature ramp conditions or at isothermal temperatures distant from ambient conditions. This test method does not address the temperature effects on mass calibration
- On Manual calibration, TA Instruments uses a zero tare, then a 100mg and 1000mg mass standards

25



25

Calibration – Weight (Auto)

Weight
Balance is steady

-0.0052 mg

Auto **Semi-Auto**

Calibration Verification

Verify automatically after calibration

Verification Criteria: Weight \pm 0.10 %

	Pan Number	Weight
Calibration Fixture 1	11	345.022 mg
Calibration Fixture 2	12	443.936 mg
Calibration Fixture 3	13	1246.285 mg

Start

Schedule


Status:
Ready

26



26

Calibration – Weight (Manual)



Weight
Balance is steady
-0.00052 mg

Auto | Semi-Auto

Calibrate Zero Range

1. Place an empty pan of the same type and size as the sample pan, on the tare side of the balance.
2. Place an empty sample pan on the loader
3. Tare the empty sample pan by clicking Tare. It may take several minutes to complete the zero range.

Tare

Calibrate Weight Range

Enter the mass of the calibration weight

1. Place a known weight in the sample pan then enter the weight in the box:
2. Click Calibrate to load the sample pan with the known weight and calibrate the weight range. It may take several minutes to complete
3. Place the Larger weight into the pan and enter the weight in the box
4. Click Calibrate to load the sample pan with the known weight and calibrate the weight range. It may take several minutes to complete

100.000 mg
Calibrate

1000.000 mg
Calibrate


27



27

Calibration – Weight (Manual)

Simply reload the mass to verify – mass difference of $\sim 0.005\%$



Weight
Balance is steady
100.1517 mg

28



28

ASTM E 1582 - Calibration of Temperature Scale for TGA

- The standard describes two methods by which the TGA can be calibrated for temperature; by melting point or magnetic transition. The most common approach for a TGA would be the magnetic transition approach
- Curie Point Temperature - that temperature where the material loses its magnetic susceptibility - defined as offset point
- Paramagnetic - a material that is susceptible to attraction by a magnet
- Temperature Calibration points are determined by comparing the measured melting onset temperature to the literature value
- TA Instruments software allows for up to 5 temperature calibration points
 - Generally, these should bracket the temperature range of interest for subsequent samples


29



29

TGA: Temperature Calibration

Vertical Balance Configuration



30




30

SDT: Temperature Calibration

Horizontal Balance Configuration

Attraction of sample to magnet, results in initial weight loss, which is lost at the Curie Transition Temperature



31

Curie Temperature Reference Materials

- International Confederation for Thermal Analysis and Calorimetry (ICTAC) developed a set of six certified and traceable Curie temperature reference materials for the calibration of TGA
 - Alumel 153°C
 - Nickel 358°C
 - Ni83Co17 555°C
 - Ni63Co37 747°C
 - Ni37Co63 931°C
 - Cobalt 1116.0°C
- The materials permit temperature calibration in about 200 °C intervals over the range of 150 to 1120 °C
- TA Instruments is the exclusive worldwide distributor for these Curie point materials

32



32

Calibration – Temperature

Calibration

Calibration Data | Temperature Calibration | DTA Signal Setup | Weight Calibration

Temperature Calibration Setup

Pan Type: Platinum (100 uL)

Operator:

Project:

Notes:

Insert Isothermal:

Ramp: °C/min

Calibration Experiments

Reference Material	Material Type	Reference Temperature	Lower Limit	Upper Limit	Pan Number
Nickel	Cure Point	358.2	258.2	458.2	1

Add Experiment

Calibration Perform Verification after Calibration

Verification Perform Calibration if Verification fails

Verification Criteria: Temperature \pm °C

Verification Experiments

Reference Material	Material Type	Reference Temperature	Lower Limit	Upper Limit	Pan Number
Nickel	Cure Point	358.2	258.2	458.2	1

Add Experiment

33



33

Verification – Temperature

Experiments

Design Run | Design View (0) | Schedule

Run 1 in Design View

Sample

Sample Name:

Pan Number:

Standard: Nickel

Pan Type: Platinum (100 uL)

Operator:

Project:

Notes:

File Name: C:\ProgramData\TA Instruments\TRIOS\Data\Default.tri

Procedure

Test: Temperature Verification

Name: Temperature Verification

Template | Segments


No.	Description
1	Equilibrate 258.2 °C
2	Electromagnet On 50.0 %
3	Ramp 20.00 °C/min to 458.2 °C

34



34

Calibration – Temperature



35



35

Calibration – Temperature

From the 'Results' Pane, Choose both the Alumel and Nickel Runs

The screenshot shows the TA Instruments software interface. The main window displays a graph of Weight (%) vs Temperature (°C) for an 'Alumel' run. The graph shows a sharp drop in weight percentage around 140°C. The interface includes several panels:

- File Manager:** Shows a list of runs including 'Alumel' and 'Nickel'.
- Control Panel:** Displays various parameters such as 'Method Time', 'Temperature', 'Sample Temperature', and 'Balance Purge'.
- Log:** Shows a table of experimental runs with columns for 'Day', 'Code', and 'Log Time'.
- General:** Contains settings for 'Total Flow', 'Balance', and 'Blending'.


36



36

Calibration – Temperature

Choose 'Calibrations' and
'TGA Temperature'




37



37

Calibration – Temperature

Dialogue Box
Shows Summary of
Calibration –
Choose 'Save and
Apply' if Data is OK




38




38

TGA – Temperature Calibration



39

TGA – Temperature Verification



40

Calibration – Temperature: Verification

No.	Description
1	Equilibrate 100.00 °C
2	Isothermal 1.0 min
3	Electromagnet On 50.0 %
4	Ramp 10.00 °C/min to 200.00 °C
5	Mark End of Cycle
6	Data Off
7	Electromagnet Off 0.0 %
8	Equilibrate 300.00 °C
9	Data On
10	Isothermal 1.0 min
11	Electromagnet On 50.0 %
12	Ramp 10.00 °C/min to 400.00 °C

Alumel


Nickel

41



41

Curie Standards with ICTAC traceability



42



42

Calcium Oxalate “Standard” Analysis




- Although Calcium Oxalate is not generally accepted as a “Standard Material,” it does have practical utility for INTRA-laboratory use
- Carefully control the experimental conditions; i.e. pan type, purge gases/flow rates, heating rate
- Particularly control the amount (~5mg) and the particle size of the sample and how you position it in the pan
- Perform multiple runs, enough to do a statistical analysis
- Analyze the weight changes and peak temperatures and establish the performance of YOU and YOUR instrument
- When performance issues come up, repeat the Calcium Oxalate analysis

43



43

Calcium Oxalate Decomposition




44



44

Calcium Oxalate Repeatability

Overlay of 8 runs, same conditions



45



45

Calcium Oxalate Repeatability

Run #	Transition 1		Transition 2		Transition 3	
	Wt Change %	Peak Temp °C	Wt Change %	Peak Temp °C	Wt Change %	Peak Temp °C
1	12.13	156.68	18.78	493.37	29.62	684.33
2	12.22	153.60	18.75	494.17	29.56	680.43
3	12.20	155.40	18.76	495.6	29.63	684.11
4	12.21	155.58	18.77	495.98	29.69	688.11
5	12.21	154.05	18.75	494.72	29.54	684.28
6	12.20	154.91	18.73	495.62	29.58	684.83
7	12.21	155.09	18.77	494.71	29.61	683.92
8	12.20	153.52	18.77	493.84	29.57	681.85
Ave	12.20	154.85	18.76	494.75	29.60	683.98
Std Dev	0.028	1.08	0.016	0.93	0.048	2.24
Theoretical	12.3		19.2		30.1	
Accuracy	0.8%		2.3%		1.7%	
Precision	0.2%		0.1%		0.2%	

46



46

SDT Calibration and Verification

- Calibrate upon initial installation
- Re-calibrate anytime the beam set, experimental heating rate, or purge gas is changed
- Types of calibration available:
 - Weight Calibration: (TGA weight signal)
 - DTA Signal Setup: Analyzing the Delta T signal data
 - Temperature (Melting point or curie point standards as in TGA. Commonly use melting point standards)
 - DSC Heat Flow
 - MDSC Reversing Heat Capacity (SDT 650)

47



47

SDT Calibration and Verification

- DTA signal:
 - Not required when using the SDT as a DSC-TGA
 - This run usually utilizes the same baseline run obtained for TGA Weight Calibration
- Heat flow and cell constant:
 - Based on analyzing the heat capacity curve for sapphire over the range 200 to 1500°C. Three experimental runs are required: two runs to generate the heat flow curve and another run to refine that calibration through cell constant calibration using a known metal standard (zinc, for example)
- MDSC Reversing heat capacity:
 - A heat capacity calibration curve is generated by running a sapphire sample over a desired temperature range using appropriate modulated conditions. The collected Reversing Heat Capacity curve is calibrated against the true value of the heat capacity of sapphire over the experimental temperature range

48



48

Instrumental Considerations



49

Instrument Hardware and Gas Selection Considerations

- Gas Delivery Module and Mass Flow Controllers
 - The gas 1 port purges both sample and balance areas
 - Gas 1 should be an inert gas (N₂, He, Ar)
 - The gas 2 port is used when a different purge gas is required, or gas switching is used
 - Typically this is air or O₂
 - Gas type is assigned to Mass Flow Controller in the Instrument section of the control software and chosen before on the setup page.

Gases Typically used on TGA/SDT

- Nitrogen - inert, inexpensive and readily available
- Helium - inert, commonly used on TGA-MS
- Argon - inert
- Air/Oxygen - used when studying oxidative stability of materials, can sometimes improve resolution of weight loss events




50



50

Blending Gas Delivery Module

- For use with TGA 550, TGA 5500 and SDT 650
- Allows blending two gases as main sample purge for a test. Nitrogen, helium, argon, oxygen, air, carbon dioxide, carbon monoxide, and forming gas (a blend of 4% hydrogen with 96% nitrogen) may be blended



51


Four gases can be used in a test



51

Blending Gas Delivery Module

- Balance and sample purge flow will depend on the instrument:



Back panel of the Blending GDM.

TGA 550	40	60	Balance/Sample
TGA 5500	25	25	
SDT 650	100	100	


52



52

Heat Exchanger - TGA


- The heat exchanger contains a liquid reservoir that supplies the instrument with coolant to dissipate heat from the furnace
- The coolant exits the heat exchanger through the supply line, circulates to the furnace, and returns to the reservoir via the return line
- Check the level and condition of the heat exchanger coolant periodically (about 3 months)
- To clean: empty old water, fill with distilled water
 - Add TA Instruments TGA Conditioner
 - (P/N 952377.901) (algae growth suppressor)
- For Q series, after filling, in software choose "Control \ Prime Exchanger"



53

53

Test for Oxygen Contamination of N2 Purge Gas




54



54

Effect of Oxygen on Copper Oxalate



55

Baseline Performance Verification

- A good way to quantify how well the TGA is working
- Especially important for measuring small weight losses associated with volatilization or small amounts of residue
- Run clean, empty, tared pan, over temperature range of interest, at desired heating rate
- Plot weight in μg vs. temperature
- Dynamic drift should be less than 10 μg for the Discovery TGA 5500, and Discovery TGA and less than 50 μg on the Discovery TGA 550/55 & Q Series TGA's when using platinum pans and 20°C/min heating rate

57

TA

57

TGA: Factors Influencing Baseline

- Stability of table
- Hang down wire / beam condition
- Hang down tube condition
- Leveling of TGA
- Cleanliness of the furnace
- Purge gas flow rates


58



58

Discovery TGA 5500 Baseline Performance

Empty Platinum Pan Baseline Scan @ 20°C/min




59



59

Effect of Purge Gas Flowrate on Baseline




60



60

Effect of Purge Gas Flowrate on Baseline




61



61

Balance Sensitivity- 120 μg Calcium Oxalate




62



62

Sample: 27 μg Sodium Tartrate




63



63


TGA: High Sensitivity Volatiles Analysis



64

TGA: Sample Preparation


- Sample mass
 - 10-20mg for most applications
 - 50-100mg for measuring volatiles or residues
- If a TGA has a baseline drift of +/-25mg then this is 0.25% of a 10mg sample



65

65

TGA particle size matters, NaCl decrepitation



66



66

TGA: Sample Preparation

- Use brass tweezers to eliminate static effects
- Tare a clean sample pan before every run
- Distribute sample evenly over bottom of pan
- Liquid samples - use hermetic pan with a pin-hole lid

67




67

TGA: Sample Pans - Types/Sizes



Deep-walled pans are good for larger mass and low-density materials




68



68

TGA: Sample Pan Selection

- Platinum (useful for most materials)
 - Easy to clean
 - Nonporous
 - Can alloy with most metals
- Alumina (Ceramic)
 - Corrosives/Inorganics
 - Large samples
 - Porous, can be easily contaminated
- Aluminum (TGA) (designed for one-time use)
 - Lower cost, disposable
 - Lower temperature limit ($\leq 600^\circ\text{C}$)



69



69

SDT: Sample Pans - Types/Sizes

- Platinum:

- 40 mL
- 110 mL


- Alumina:

- 40 mL
- 90 mL



70

Sealed Aluminum Pans




71



71

Sealed Aluminum Pans and Punching




72



72

TGA: Sample Pan Cleaning

- All sample pans are reusable (except Aluminum)
- If using platinum or alumina pans, a flame torch can be used to burn off organic residue. (do not flame Aluminum pans)
- Scrape off remaining ash (DSC fiberglass brush)
- Swab out with an organic solvent such as acetone or alcohol. Let it dry out before using it



73



73

TGA: Environment Considerations

- Avoid areas near heater or air conditioner ducts
- Avoid tables with drawers or those near a door
- For optimum results, use a marble table

74



74

Experimental Methods




75

75

Experiment Setup

- Experiments are created in the Running Queue or the Design View
- They are launched from the Running Queue!



76



76

Experiment Setup

- Sample information is entered here

Sample

Sample Name

Pan No.

Pan Type

Operator

Project

Notes



File Name: C:\Documents and Settings\All Users\Application Data\TA Instruments\TRIOS\Data\Default.tri


77





77

Experiment Setup

Procedure  

Test 

Name

Template  Segments 

Heating Rate °C/min

Final Temperature °C


Switch to gas 2 at °C

78



78

Experiment Setup

Advanced 

Beginning of Test:

Start Experiment After Weight Stabilization

Uncheck weight stabilization for volatile materials

End of Test:

Enable Air Cool

Uncheck for EGA to prevent background gases

Air Cool Until Temperature Is Below °C

End of Test Delay minutes


79



79

Segment Statements

- The logic of the instrument control software is based upon segment statements which the user enters during the design of the experiment
- These segments are executed during the experiment.




80



80

Method Segments




81



81

Method Design: TGA Segment List

- The Ramp segment heats or cools the sample at a fixed rate until it reaches the specified temperature, producing a linear plot of temperature versus time
- The Equilibrate segment heats or cools the furnace to the defined temperature, stabilizes the furnace at that temperature, then continues to the next segment
- The Select Gas segment controls the switching of gas between Gas 1 and Gas 2 for an instrument with a gas delivery module. This segment is used to synchronize gas switching at a specific time or temperature in an experiment






82



82

Segments and Descriptions






Segment	Description
 Abort	<p>The Abort segment skips over the next segment when specified limit conditions are met.</p> <ul style="list-style-type: none"> • If the limit is reached at the beginning of a segment, then that segment is skipped and method execution continues with the next segment. • If the limit is reached during the execution of a segment, then the remaining portion of the segment is skipped. <p>• NOTE: The Abort segment is generally followed by a Ramp or Isothermal segment.</p> <p>Example (DSC):</p> <ol style="list-style-type: none"> 1. Equilibrate at 200°C 2. Abort next segment if mW>1 3. Isothermal for 100 min
 Balance Flow Rate (Discovery TGA only)	<p>This segment is used to alter the rate of flow of the selected gas to the balance.</p> <p>Example: Flow rate 50 mL/min</p>
 Blend Gas	<p><i>Applicable to Blending GDM instruments only:</i> The Blend Gas segment allows you to select the input gas for Channel A (Gas 1 or Gas 2), a percentage to blend with Channel B (Gas 3 or Gas 4), and which input Gas to use for Channel B.</p> <ul style="list-style-type: none"> • NOTE: Minimum controllable flow rate is 5 mL/min. Take this into account when specifying percentage. It may be necessary to increase overall sample flow. <p>Example: Blend Gas 1 at 60% with Gas 3</p>

83



83

Segments and Descriptions







Data 	<p>The Data segment controls data collection during the experiment. If a Data segment is not used, data storage is automatically initiated by the first Ramp, Isothermal, or Step segment that appears in the method.</p> <p>Example: Data Storage: On</p>
Electromagnet (Discovery TGA only) 	<p>The Discovery TGA has a magnetic coil surrounding the furnace. The Electromagnet segment allows you to apply a magnetic field during an experiment so that temperature calibration using Curie point standards may be performed.</p> <p>Example: Electromagnet: On Ramp 10°C/min to 250°C</p>
Equilibrate 	<p>The Equilibrate segment heats or cools the furnace to the defined temperature, stabilizes the furnace at that temperature, then continues to the next segment. This segment does not automatically start data collection.</p> <p>Example: Equilibrate at 200°C</p>
Event 1 / Event 2  	<p>The Event segment controls the external event relay through the event jack on the back of the instrument. This is used to synchronize control of additional hardware through the method.</p> <p>Example: Event 1: On</p>

84



84

Segments and Descriptions

 Initial Temperature	<p>The Initial Temperature segment heats or cools the furnace to the defined temperature, stabilizes the furnace at that temperature, then holds the temperature until the experiment is continued by clicking OK on the TRIOS dialog box, or by selecting Start on the instrument display or instrument keypad. This segment does not automatically start data collection.</p> <p>Example: Initial Temperature 200°C</p>
 Isothermal	<p>The Isothermal segment holds the sample at the current temperature (as programmed by the previous segment) for a defined period of time. This segment automatically turns on data collection, except when preceded by a Data OFF segment.</p> <p>Example: Isothermal for 10 min</p>
 Jump	<p>The Jump segment instantly changes the set point temperature, causing ballistic changes in the sample temperature. This segment then allows the immediate execution of the next segment (which is usually the Isothermal segment). Note that large temperature overshoots may result from the use of this segment. This segment does not automatically start data collection.</p> <p>Example: Jump to 200°C</p>
 Mark End	<p>The Mark End segment places a marker in the data for use by the data analysis programs. In general, markers provide quick parsing of data to separate experimental segments (i.e., the heat-cool cycle).</p> <p>Example: Mark end of cycle 0</p>
 Mass Flow	<p>The Mass Flow segment alters the rate of flow of the selected gas when an instrument is equipped with a Gas Delivery Module (GDM).</p> <p>Example: Mass Flow 50 mL/min</p>
 Modulate Temperature	<p><i>Available for Modulated Instruments Only:</i> This segment allows you to enter the modulation temperature amplitude and period (frequency) parameters that will be used with subsequent ramp or isothermal segments.</p> <p>Example:</p>

85



85

Typical Methods

- Ramp (heating) experiment:
 - Ramp 20°C/min. to 800°C
- Ramp and switch gas (carbon black content, residue)
 - Ramp 20°C/min. to 650°C
 - Select gas: 2
 - Ramp 20°C/min. to 1000°C

86



86

What if I need help?

- TA Tech Tips
 - <http://www.youtube.com/tatechtips>
- TA Instruments Applications Helpline available from the TA website
 - <http://www.tainstruments.com/support/applications/applications-hotline/>
- Check out our Website
 - <http://www.tainstruments.com/>

87



87

TA Instruments –Waters LLC
www.tainstruments.com

