

# Discovery DSC 25, DSC 250, DSC 2500



## Site Preparation Guide

# Table of Contents

|  |       |
|--|-------|
| Table of Contents .....                | 2     |
| Ideal Setup .....                      | 3     |
| System Components.....                 | 4     |
| Instrument Measurements .....          | 5     |
| Utility Requirements.....              | 6     |
| Computer Requirements.....             | 7     |
| Accessories .....                      | 8     |
| Refrigerated Cooling System (RCS)..... | 8–9   |
| Finned Air Cooling System (FACS) ..... | 10    |
| Liquid Nitrogen Pump (LN Pump).....    | 11-12 |
| Photocalorimeter Accessory (PCA).....  | 13    |
| Site Preparation Checklist .....       | 14    |
| TA Instrument Offices.....             | 15    |

# Ideal Setup



## IDEAL PLACEMENT AND BENCH MEASUREMENTS

Select a location with a rigid laboratory bench that is level and is in a vibration-free environment.



Bench length: 183 cm (72 in)

Bench depth: 76 cm (30 in)

Distance from the wall: 30.5 cm (12 in) min.




The casters on the LN Pump dewar are 61 cm (24 in) x 61 cm (24 in). Allow for 1–3 ft of space between the DSC and the LN Pump depending on how the supply/return line is oriented.

# System Components



## MAIN SYSTEM COMPONENTS



- A. Cooling Accessory (RCS shown)
- B. Instrument
- C. Computer (Controller)

# Instrument Measurements



## DSC WITH AUTOSAMPLER




Height: 61 cm (24 in)

Width: 53 cm (21 in)

Depth: 51 cm (20 in)



## DSC WITHOUT AUTOSAMPLER



Height: 43 cm (17 in)

Width: 53 cm (21 in)

Depth: 51 cm (20 in)



Weight\*: 22 kg (48 lbs)

\*Includes Autosampler, Autolid, and FACS

# Utility Requirements



## POWER

- 100–240 VAC, 50/60 Hz, 600 W
- Safety ground per local regulation



## GAS

Cell purge gases: air, nitrogen, oxygen, argon, helium

| Description                                      | Requirement  |
|--|--|
| Cell/base purge gas pressure:                    | 100–140 kPa (10–20 psig)   |
| Cooling gas (air) pressure for <b>FACS</b> :     | 170 kPa (25 psig max)  |
| Cooling gas (nitrogen) pressure for <b>RCS</b> : | 170 kPa (25 psig max)  |
| Conditions                                       | <ul style="list-style-type: none"><li>• Must be dry</li><li>• Must be free from oil and dirt</li></ul> |
| Cooler   | Use dry nitrogen as the base purge gas when using a cooler.  |
| Other  | 1/8" polyethylene tubing and fittings are supplied in the accessory kit                                |



RCS90

# Computer Requirements



## HARDWARE REQUIREMENTS

| Description       | Requirement  |
|-------------------|--|
| Processor         | <ul style="list-style-type: none"><li>• Intel® Core™ 2 Duo or better</li><li>• 2.93 GHz with 3 MB L2 cache</li></ul>   |
| Memory            | ≥ 16 GB RAM  |
| Hard drive        | ≥ 80 GB free space <ul style="list-style-type: none"><li>• 700 MB required for Full version of TRIOS</li><li>• 400 MB required for Lite version of TRIOS (without Online help)</li></ul> |
| DVD               | ≥ 48x CD-ROM or DVD  |
| Screen resolution | Required: 1280 x 1024 with 24-bit colors<br>Recommended: 1920 x 1080 with 24-bit colors  |
| Graphic memory    | 128 MB   |
| Screen (LCD) size | Required: 19" or greater<br>Recommended: 24" wide screen   |

# Computer Requirements



## SOFTWARE REQUIREMENTS

| Item             | Requirement   |
|------------------|---|
| Operating System | <ul style="list-style-type: none"><li>• Windows 7, 8, 10 Ultimate, Enterprise, &amp; Professional</li><li>• Home version not supported</li></ul>                |
|                  | 64-bit version for TRIOS v4.5 and up  |
| Other            | Microsoft Operating System Service Pack   |
| Other            | Must be up-to-date  |
| Browser          | Internet Explorer   |
| Network          | <i>TA Instruments is not responsible for resolving issues associated with connections to your corporate network.</i>  |
| Conflicts        | <i>TA Instruments is not responsible for resolving hardware/software conflicts created by the addition of third party hardware or software to the computer.</i> |



# Accessories



## REFRIGERATED COOLING SYSTEM MEASUREMENTS











|        | RCS120            | RCS90             | RCS40            |
|--------|-------------------|-------------------|------------------|
| Height | 88 cm (34.6 in)   | 46 cm (18 in)     | 26 cm (10 in)    |
| Width  | 35.6 cm (14.5 in) | 26 cm (10 in)     | 26 cm (10 in)    |
| Depth  | 56 cm (22 in)     | 51 cm (20 in)     | 51 cm (20 in)    |
| Weight | 102 kg (225 lbs)  | 47.7 kg (105 lbs) | 24.8 kg (55 lbs) |

# Accessories



## REFRIGERATED COOLING SYSTEM REQUIREMENTS

### Requirements

|   | RCS120   | RCS90                                   | RCS40  |
|---|--|---|--|
|   | 230 VAC/8.5 A/50 Hz<br>230 VAC/7.5 A/60 Hz   | 120 VAC/12 A/60 Hz<br>220 VAC/6 A/50 Hz | 120 VAC/6.25 A/60 Hz<br>220 VAC/4 A/50 Hz  |
|    | Uses NEMA 5-15 plug  | Uses NEMA 5-20P plug                    | Uses NEMA 5-15 plug  |
|   | <br>NEMA 5-20P  |   | <br>NEMA 5-15 |
|   | <ul style="list-style-type: none"> <li>RCS90/RCS40: Place the RCS90 and RCS40 on a table separate from the laboratory bench. If no table is available, place the RCS on the bench to the left of the instrument.</li> <li>RCS120: RCS120 must be kept on the floor</li> </ul>  |   |  |
|  | <ul style="list-style-type: none"> <li>A base and cooling purge (nitrogen) is required in addition to the standard cell purge</li> <li>Use 99.999% pure nitrogen or LN boil-off gas to reduce moisture</li> <li>New or recently used calibrated regulator is recommended</li> <li>Make sure tubing is cut cleanly and squarely on the ends. Use of the Legris Tubing Cutter #3000-71-00 is recommended </li> <li>Leak check all tubing</li> </ul> |   |  |
|   |  Do not use Tygon® tubing due to its high moisture permeability   |   |  |
|  | Customer-supplied: <ul style="list-style-type: none"> <li>Regulator</li> <li>Moisture trap (P/N 200266.001) to prevent moisture build-up</li> </ul>  |   |  |



Circulator



Power



Cooling



Gas



LN<sub>2</sub>



Fluid



Light



Hardware



Software



Temp



Lab



Customer

# Accessories



## FINNED AIR COOLING SYSTEM REQUIREMENTS

### Requirements



- Cooling gas (air) maximum air pressure: 25 psig (170 kPa gauge)
- Use standard grade nitrogen and clean house air
- Leak check all tubing
- Make sure tubing is cut cleanly and squarely on the ends. Use of the Legris Tubing Cutter #3000-71-00 is recommended



#### Recommendations:

- Use a filter
- Use a new or recently calibrated regulator



Customer-supplied: Regulator



Circulator



Power



Cooling



Gas



LN<sub>2</sub>



Fluid



Light



Hardware



Software



Temp



Lab



Customer

# Accessories



## LN PUMP MEASUREMENTS



Height: 122 cm (48 in)

Width: 86 cm (34 in)

Depth: 86 cm (34 in)

Weight EMPTY: 50 kg (110 lbs)

Weight FULL: 92.5 kg (204 lbs)

# Accessories



## LN PUMP REQUIREMENTS

### Requirements



24VDC using universal power supply (refer to the Discovery LN Pump Getting Started Guide)

#### Requirements

- Cooling gas (nitrogen or LN boil-off) maximum pressure for use with the LN2P = 170 kPa gauge (25 psig)
- Low pressure Liquid Nitrogen dewar required




Do not use Tygon® due to its high moisture permeability

- Use new or recently serviced/calibrated regulator
- Use 99.999% pure helium to reduce moisture build-up in the cell



#### Recommendations

- Helium gas recommended for cell purge via the GAS 2 port
- Make sure tubing is cut cleanly and squarely on the ends. Use of the Legris Tubing Cutter #3000-71-00 is recommended 
- Leak check all tubing
- Use the gas dryer (P/N 200266.001) to pre-dry and indicate unsatisfactory moisture levels
- Use the purge gas purifier (P/N 970425.901) to achieve a dew point of -180°C



Circulator



Power



Cooling



Gas



LN<sub>2</sub>



Fluid



Light



Hardware



Software



Temp



Lab



Customer

# Site Preparation Checklist



## PHOTOCALORIMETER ACCESSORY MEASUREMENTS



The PCA cannot be used with the DSC 25.



Height: 15 cm (6 in)

Width: 28 cm (11 in)

Depth: 44 cm (17 in)

Weight: 9.4 kg (21 lbs)



## PHOTOCALORIMETER ACCESSORY REQUIREMENTS

### Requirements



Same general requirements as DSC. See pages 5–7.



Circulator



Power



Cooling



Gas



LN<sub>2</sub>



Fluid



Light



Hardware



Software



Temp



Lab






Customer

# Site Preparation Checklist



## Discovery DSC 25, DSC 250, DSC 2500

|   |   |
|---|---|
|    | <p>Sufficient bench space for instrument, computer, and cooling accessory (if needed)</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Length: 183 cm (72 in)</li> <li><input type="checkbox"/> Depth: 76 cm (30 in)</li> </ul>  |
|    | <p>Instrument power is 100–240 VAC, 50/60 Hz, 600 W</p>   |
|    | <p>Purge gas– Air, nitrogen, oxygen, argon, or helium</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Base purge gas pressure is 100–140 kPa (10–20 psig)</li> <li><input type="checkbox"/> Cooling gas (air) pressure for FACS is 170 kPa (25 psig max)</li> </ul> <p style="text-align: center;">OR</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Cooling gas (nitrogen) for RCS or LN Pump is 170 kPa (25 psig max)</li> <li><input type="checkbox"/> Use 99.999% pure nitrogen or LN boiloff gas to reduce moisture</li> <li><input type="checkbox"/> Moisture trap (P/N 200266.001) to prevent moisture buildup</li> </ul> |
| <p>I hereby acknowledge that all utility requirements have been met per the checklist above and that they will be ready at the agreed time of installation.</p> |   |
| <p><b><u>If all utility requirements are not met at the agreed time of installation, additional charges may be incurred for a return Service trip.</u></b></p>  |   |
| <p>_____</p> <p><i>Customer</i></p>   | <p style="text-align: center;">____/____/____</p> <p style="text-align: center;"><i>DD MM YYYY</i></p>  |
| <p>_____</p> <p><i>Company</i></p>  | <p style="text-align: center;">_____</p> <p style="text-align: center;"><i>City</i></p>   |
| <p>_____</p> <p><i>Company</i></p>  | <p style="text-align: center;">_____</p> <p style="text-align: center;"><i>State</i></p>  |
| <p>_____</p> <p><i>Company</i></p>  | <p style="text-align: center;">_____</p> <p style="text-align: center;"><i>Country</i></p>  |
| <p>Please send a signed copy of the completed checklist to your local Service representative.</p>   |   |

# TA Instruments Offices

For information on our latest products, contact information, and more, see our website at:  
<http://www.tainstruments.com>.

To find your local TA Instruments office and contact information, visit  
<http://www.tainstruments.com/contact/ta-directory/>

TA Instruments – Waters LLC  
Corporate Headquarters  
159 Lukens Drive  
New Castle, DE 19720  
USA

Telephone: 302-427-4000  
Fax: 302-427-4001  
Email: [info@tainstruments.com](mailto:info@tainstruments.com)