

# Simultaneous DSC-synchronous XRD for understanding solid-state phase transformations

Professor Simon Gaisford

[s.gaisford@ucl.ac.uk](mailto:s.gaisford@ucl.ac.uk)



## Scope

- Background to physical forms, especially in relation to pharmaceuticals
- Physical form identification by DSC
- Methods that might aid interpretation
- Hyphenated DSC techniques
- DSC-XRD
- Summary

## DSC for physical form characterisation

- Methodology simple
- Put sample in a pan, heat it then look for peaks!
- Tricky part is interpretation of the data
  
- No molecular information and different processes can give similar peaks
- How to interpret?
  - LOOK at whether events are exo- or endothermic


What's going on here?




# DSC for physical form characterisation

- Methodology simple
- Put sample in a pan, heat it then look for peaks!
- Tricky part is interpretation of the data
  
- No molecular information and different processes can give similar peaks
- How to interpret?
  - LOOK at whether events are exo- or endothermic
  - Much can be gained from altering DSC method, *especially* by reheating sample and by changing heating rate

# Reheating sample



# A metastable polymorph at two heating rates



## Additional analyses

- Can help interpretation with additional analyses
  - TGA
  - Hot-stage microscopy
  - Evolved gas analysis
  - Raman







Telford et al (2016)  
 10.1039/c6cc05006a

Paracetamol


(a)



(b)



- Data for glutaric acid
- The question is, to what events should the peaks be assigned?




## DSC for physical form characterisation

- ‘Gold standard’ analysis is single-crystal or powder XRD
- Takes a long time to record diffraction pattern
- We wanted to confirm physical form with PXRD in real-time, in the DSC pan
- Best done with a synchronous X-ray source







NODE 804


1

6

USER CONTROL

PLEASE CONTACT THE OPERATOR

UCL  
UCL  
UCL









## DSC for physical form characterisation

- Problem comes in *absolute* confirmation of form
- Can use *in-situ* Raman or NIR, but usually single-crystal or powder XRD is required
- Takes a long time to record diffraction pattern
- We wanted to confirm physical form with PXRD in real-time, in the DSC pan
- Best done with a synchronous X-ray source
  
- We were able to record a powder pattern in 4s, with a 2s pause, so 1 pattern every 1 °C

## How to interpret XRD data?







Glutaric acid:

Two enantiotropic polymorphs


Solid-solid conversion  
From  $\beta$  to  $\alpha$



By selecting specific XRD reflections, can look at loss of  $\beta$  and creation of  $\alpha$




## Sulphathiazole



Sulphathiazole:


Three enantiotropic polymorphs


Solid-solid conversion from a mixture of forms III and IV to form I





Again, can see loss of Forms III and IV and creation of Form I






(a)


Telford et al (2016)  
10.1039/c6cc05006a

# Paracetamol



# ChemComm

Chemical Communications  
[www.rsc.org/chemcomm](http://www.rsc.org/chemcomm)





ISSN 1359-7345






COMMUNICATION  
 Ian J. Scowen *et al.*  
 Stabilisation of metastable polymorphs: the case of paracetamol form III

**175** YEARS



# 10,11-dihydrocarbamazepine









IV I L





II & III I L





## Olanzapine





## Summary

- DSC excellent for seeing phase transitions, but absolute assignment to physical forms tricky
- Simultaneous XRD capable of following change in form in real-time
- TA DSC is very easy to modify for this use
- Initial experiments have shown that solid-solid transitions occur without an intermediate liquid phase
- If you have systems of interest, contact us and we will schedule them on the Diamond Beam

## Thanks to...

- Gareth Williams
- Asma Buanz
- Alex Clout
- Sean Askin
- Richard Telford and Ian Scowen
- Sally Price, Derek Tocher and Jeremy Cockcroft
- Jas Mahey and TA Instruments



**EPSRC**

Engineering and Physical Sciences  
Research Council



1949 version

2016 version



 **diamond** is a DSC's best friend